# taxon-utils之计算给定分类的谱系: lineage

## 一、taxon-utils lineage介绍

#### 功能描述:

taxon-utils lineage 计算给定分类的谱系.

#### 命令行接口:

```
$ taxon-utils lineage

Usage: taxon-utils lineage [options] <taxon.map> <tab|taxonIds>

Options:
   -c INT column for taxon translate, default: [2]
   -n print lineage for specified node, using ',' for multi taxonId.
```

#### 可选参数:

```
1 -c 整型 指定转化的列,默认为第二列;
2 -n 输出指定节点的世系
```

### 二、使用场景实例及其用法

#### 示例演示:

示例文件在: data 目录

示例文件: classify.txt.gz, taxon.map.gz

```
1 | $ zcat classify.txt.gz | head -n 6
```

```
C A01050:204:HF7FGDSXY:4:1101:12943:1016 25457 150|149 0:22 25457:4 0:7 25458:5 0:78 |:| 3:4 3343:5 8109:2 25457:5 25460:1 25457:6 25693:5 0:48 25458:2 0:1 3343:2 14694:7 3343:10 3:5 0:12 C A01050:204:HF7FGDSXY:4:1101:12337:1031 18331 150|149 0:60 31294:5 0:25 18331:3 0:23 |:| 0:93 18331:3 0:19 C A01050:204:HF7FGDSXY:4:1101:16866:1047 8364 150|150 0:17 20367:2 35825:1 0:8 3:4 20616:5 3:1 4523:5 3:3 0:9 3:5 0:19 27173:2 0:12 27173:2 0:21 |:| 0:2 2878:3 24660:5 28576:3 30576:2 3:13 0:9 2483:17 27170:3 2483:2 27170:5 3:6 8364:3 3:7 8364:3 27172:2 8364:5 27172:6 171:5 0:15 C A01050:204:HF7FGDSXY:4:1101:22742:1047 23158 150|150 0:82 23158:4 0:30 |:| 0:116 C A01050:204:HF7FGDSXY:4:1101:28664:1063 559 150|150 0:22 559:3 0:7 171:5 0:79 |:| 0:116 C A01050:204:HF7FGDSXY:4:1101:20473:1094 25986 150|150 0:45 25986:3 0:68 |:| 0:116
```

注意事项: classify 文件为 Kraken2 分类的结果,第一列为标识符: C: 可以分类的序列, U: 不能分类的序列,第三列为分类的 Taxonomy ID,可以为 NCBI 分类号或者使用 GTDB 的自定义分类号,正常我们只需要第二列和第三列。

```
1 | $ zcat taxon.map.gz | head -n 6
```

```
1
1
  1
                no rank root
                              root
2
  2
         1
                 superkingdom Archaea root
3
  3
         1
                superkingdom Bacteria
                                           root
4
         3
                phylum 4572-55 Bacteria
  4
                 phylum AABM5-125-24 Bacteria
5
  5
         3
6 6
                 phylum AB1-6 Bacteria
         3
```

#### 运行命令:

```
1 | $ zcat classify.txt.gz | grep -P "^C" | cut -f2,3 | head -n 6
```

```
1 A01050:204:HF7FGDSXY:4:1101:12943:1016 25457
  k:Bacteria,p:Eremiobacterota,c:Eremiobacteria,o:UBP12,f:UBA5184,g:Palsa-
   1478,s:Palsa-1478 sp003140215 3,56,265,1390,3343,8109,25457
2 A01050:204:HF7FGDSXY:4:1101:12337:1031 18331
  k:Bacteria,p:Gemmatimonadota,c:Gemmatimonadetes,o:Gemmatimonadales,f:Gemmatim
   onadaceae,g:Fen-1231,s:Fen-1231 sp003171215
   3,78,283,806,2237,5826,18331
  A01050:204:HF7FGDSXY:4:1101:16866:1047 8364
  k:Bacteria,p:Actinobacteriota,c:Actinobacteria,o:Actinomycetales,f:Micrococca
  ceae,g:Pseudarthrobacter 3,9,171,559,2483,8364
  A01050:204:HF7FGDSXY:4:1101:22742:1047 23158
  k:Bacteria,p:Proteobacteria,c:Alphaproteobacteria,o:Rhizobiales,f:Beijerincki
   aceae,g:Methylobacterium,s:Methylobacterium sp003173775
  3,110,175,1021,1696,7507,23158
5 A01050:204:HF7FGDSXY:4:1101:28664:1063 559
  k:Bacteria,p:Actinobacteriota,c:Actinobacteria,o:Actinomycetales
  3,9,171,559
6 A01050:204:HF7FGDSXY:4:1101:20473:1094 25986
   k:Bacteria,p:Actinobacteriota,c:Actinobacteria,o:Actinomycetales,f:Dermatophi
   laceae,g:Pedococcus,s:Pedococcus sp001426245
   3,9,171,559,1919,8203,25986
```

将第二例翻译成物种世系格式:

```
$ zcat classify.txt.gz | grep -P "^C" | cut -f2,3 | head -n 6 | taxon-utils
   lineage -c 2 taxon.map.gz -
2
  A01050:204:HF7FGDSXY:4:1101:12943:1016 25457
   root, Bacteria, Eremiobacterota, Eremiobacteria, UBP12, UBA5184, Palsa-1478, Palsa-
   1478 sp003140215
                       1,3,56,265,1390,3343,8109,25457
  A01050:204:HF7FGDSXY:4:1101:12337:1031 18331
   root, Bacteria, Gemmatimonadota, Gemmatimonadetes, Gemmatimonadales, Gemmatimonada
   ceae, Fen-1231, Fen-1231 sp003171215 1,3,78,283,806,2237,5826,18331
  A01050:204:HF7FGDSXY:4:1101:16866:1047 8364
   root, Bacteria, Actinobacteriota, Actinobacteria, Actinomycetales, Micrococcaceae,
   Pseudarthrobacter 1,3,9,171,559,2483,8364
  A01050:204:HF7FGDSXY:4:1101:22742:1047 23158
   root,Bacteria,Proteobacteria,Alphaproteobacteria,Rhizobiales,Beijerinckiaceae
   ,Methylobacterium,Methylobacterium sp003173775
   1,3,110,175,1021,1696,7507,23158
  A01050:204:HF7FGDSXY:4:1101:28664:1063 559
   root, Bacteria, Actinobacteriota, Actinobacteria, Actinomycetales 1,3,9,171,559
8 A01050:204:HF7FGDSXY:4:1101:20473:1094 25986
   root, Bacteria, Actinobacteriota, Actinobacteria, Actinomycetales, Dermatophilacea
   e,Pedococcus,Pedococcus sp001426245 1,3,9,171,559,1919,8203,25986
```

多输出两列,一列为主要分类水平: 第一列为完整的世系分类,一列为对应的数值编号,比如: root,Bacteria,Actinobacteriota,Actinobacteria,Actinomycetales,Dermatophilaceae,Pedo coccus,Pedococcus sp001426245,分类之间使用 ,分割,相对于 translate 具有完整的分类路 径。

命令行参数 1: -n 指定需要翻译的数字编号, 比如: 25986。

```
1 | $ taxon-utils lineage -n taxon.map.gz 25986
```

1 root,Bacteria,Actinobacteriota,Actinobacteria,Actinomycetales,Dermatophilacea e,Pedococcus,Pedococcus sp001426245 1,3,9,171,559,1919,8203,25986

本文材料为 BASE (Biostack Applied bioinformatic SEies ) 课程 Linux Command Line Tools for Life Scientists 材料,版权归上海逻捷信息科技有限公司所有。

Last Update: 2020-08-10 11:56 AM